

SOLARIS™
BIOTECH

By Donaldson



10
CATALOGUE





10

single & parallel fermenter bioreactor

IO is a fully electrical mini fermenter/bioreactor with a fast and accurate thermoregulation system (without water circulation).

“Having the opportunity to extract real time process data with the lowest step time of 1 second remotely offers to our customer a powerful monitoring tool.”

APPLICATIONS

Our products are utilized across various sectors including **pharmaceuticals, food production, agriculture, biomanufacturing, bioplastics, biofuels, cosmeceuticals**, catering to both research and development and lab production needs

OUR PRODUCT	4
PRODUCT DESCRIPTION	6
SOFTWARE LEONARDO	8
TECHNICAL DATA	10
GET IN TOUCH	12

Our Product IO

This catalogue describes Solaris IO.
For supervisory control and data acquisition
Leonardo software is included.



THE SYSTEM

The system consists of a benchtop, pre-assembled fermenter/bioreactor, supplied with all necessary tubes, valves and instruments, automation, control panel (HMI).

THE DESIGN

The system is designed for aerobic and anaerobic cultivations/fermentations, closed aseptic operations. The control is based on a SCADA control system.

Customizable Configuration

IO is completely electrical.

The thermoregulation (both heating and cooling) is performed through a Peltier Cell, placed on the bottom of the fermenter/bioreactor.

This avoids water circulation (no water source is needed in the lab).

Multiple operations: up to 24 parallel units.



BASIC
RESEARCH




EDUCATION



PROCESS
DEVELOPMENT &
OPTIMIZATION



SCALE UP &
SCALE-DOWN
STUDIES



SMALL
PRODUCTION



Fully electric:
no water circulation.

Up to 24 units managed with one HMI with innovative PARALLEL process control.

Single-wall borosilicate glass vessel, with thermoregulation performed through a Peltier Cell.

Different configurations available for microbial and cell culture applications, with the choice of Rushton/Marine/Pitched-Blade impellers and fluted/L-shaped sparger.

Modbus digital sensors reduce background noise and guarantee quick response time.

Suitable for batch, fed-batch and continuous processes.



Product Description



Different gas mixing strategies with up to 5 TMFC.

Powerful and accurate (1 RPM) brushless motor.

Wide range of measurement and control options.

Optional integration of up to 4 analog input/output connections, choosing between 0-10 V and 0-20 mA/4-20 mA (e.g. pumps or valves with power supply independent from Solaris electrical cabinet).

Extremely compact system maximizes lab space.

Additional parameter in modular external boxes for future PCS upgrade including dCO_2 , cell density, weight, peristaltic pumps.


Leonardo Software



The innovative SCADA software Leonardo has a smart and user-friendly controller designed to provide a high level of automated management of the fermentation/cultivation processes.


Full version included with the equipment. Up to 24 units managed in parallel with a unique HMI (24"). Data extraction in .csv format.

Remote control and access via PC, tablet or smartphone, with QR code scanning or via dedicated portal.



SYNOPTIC

- Real time 3D view
- Parallel control
- Manual control



LOGIC PARSER

- Customized logic functions
- Parallel logic blocks and functions



REMOTE CONTROL

- Unlimited number of profiles editor
- Unlimited number of devices to be associated



TRENDS

- Custom acquisition time
- Up to 6 values simultaneously display
- Automatic graph comparison



WORKFLOW

- Custom phase manager
- Parallel visualization
- Cascade settings
- Peristaltic pumps function assignable from software



CALIBRATION

- Up to three-point calibration
- Simultaneous calibration values for parallel work

Technical Data

VESSEL		
Solaris Code	IO 200	IO 1000
Total volume (mL)	200	1000
Ratio H/D	1.5/1	2.5/1
Min working volume (mL)	70	360
Max working volume (mL)	150	700
Max temperature	56 °C	
Operating pressure	< 0.5 bar(g)	
Headplate ports	3 ports PG13.5 (sensors, condenser, multi feed, or plugs) 2 ports M11 (sparger, harvest, LEDA, or plugs) 3 ports M12 (anti-foam probe, single feed/overlay, gas out, or plugs)	5 ports PG13.5 (sensors, condenser, adjustable level, multi feed or plugs) 5 M12 (sparger, harvest, gas out, LEDA, antifoam probe, fixed level probe, single feed/overlay or plugs)
Materials	Borosilicate Glass and AISI 316L	
SENSORS LENGTH (mm)		
Length	120	225
DIMENSIONS FOR AUTOCLAVE (with Condenser)		
Height (mm)	320	420
Diameter (mm)	170	170
STIRRINGS		
Drive and speed (rpm)	Brushless Motor, 1 - 2000 rpm	
Impellers	Select from: Rushton, Marine, Pitched blade	
THERMOREGULATION		
Control	PID Control - Accuracy 0.1 °C - Peltier cell	

GAS CONTROL & GAS MIXING

Sparger and overlay gas control	TMFC
Gas mixing (Air, CO ₂ , O ₂ , N ₂)	1 TMFC (included in entry level) + up to 4 additional TMFC
Sparger type	Select from: Flute type or L-shaped
Gas out	0.2 µm filter

PERISTALTIC PUMPS

Quantity & type	2 Watson Marlow type 400 F/A 2 Watson Marlow type 114 FD/DV
-----------------	----------------------------------------------------------------

CONTROLLER

Master control module (w x d x h)	From 1 to 24 units - 35 cm x 35 cm x 35 cm
HMI with Leonardo software	Touch screen PC, 24" color monitor

pH

Sensor	Digital sensor
Sensitivity	57 to 59 mV/pH
Control system	Measuring resident in Leonardo software
Control range	0 - 14
Operating temperature	Up to 130 °C
Pressure range	0 - 6 bar(g)
Actuator	Cascade to peristaltic pumps for the addition of acid/base solutions or gas (CO ₂)

dO₂

Sensor	Digital Optical sensor
Accuracy	1 ± 0.05 %-vol, 21 ± 0.2 %-vol, 50 ± 0.5% -vol
Control system	Measuring resident in Leonardo software
Control range	0 - 300 % air saturation
Operating temperature	Up to 130 °C
Pressure range	0 - 12 bar(g)
Actuator	Cascade to RPM, gas control, feedings

REDOX (ORP)

Sensor	Digital sensor
Control system	Measuring resident in Leonardo software
Control range	± 1500 mV
Operating temperature	Up to 130 °C
Pressure range	0 - 6 bar(g)

TECHNICAL DATA

ANTIFOAM/LEVEL

Sensor	Conductivity based
Control	Measuring resident in Leonardo software

CONDUCTIVITY

Sensor	Digital sensor
Accuracy	$\pm 3\%$ at 1 $\mu\text{S}/\text{cm}$ to 100 mS/cm , $\pm 5\%$ at 100 to 300 mS/cm
Control system	Measuring resident in Leonardo software
Control range	1 - 300 000 $\mu\text{S}/\text{cm}$
Operating temperature	Up to 130 °C
Pressure range	0 - 20 bar(g)

dCO₂

Sensor	Digital sensor
Accuracy	$\pm (10\%$ of the reading + 10 mbar)
Control system	Measuring resident in Leonardo software
Control range	0 - 200 % saturation
Operating temperature	Up to 130 °C
Pressure range	0 - 4 bar(g)

CELL DENSITY

Sensor	Digital sensor
Control system	Measuring resident in Leonardo software
Operating temperature	Up to 140 °C
Pressure range	Up to 10 bar(g) (150 psi)
TCD measuring range	0 - 4 AU
VCD measuring range	0.0 to 400 pF/cm

WEIGHT

Sensor	Digital balance
Accuracy	$\pm 0.01\text{ g}$ (3.6 kg) , $\pm 0.1\text{ g}$ (36.0 kg)
Control system	Measuring resident in Leonardo software

PERISTALTIC PUMPS

WM 120 U Brushless	1 - 100 rpm
--------------------	-------------

Get in touch



Solaris Biotech By Donaldson



info@solarisbiotech.com

**SALES &
PRODUCTION**

Via Bachelet, 58/89
46047
Porto Mantovano
Mantova - Italy

**OUR OFFICES
SHOWROOMS**





www.solarisbiotech.com

Your official Solaris Biotech distributor:

KREIENBAUM Neuroscience GmbH
Robert-Koch-Straße 9
40764 Langenfeld
Tel.: +49 (0) 2173 39927-0
E-Mail: info@kreienbaum-neo.de

