CPAPNose

KREIENBAUM NEOSCIENCE

DE Nasenschutz – Bedienungsanleitung
GB Nose protection – Instructions for Use
FR Protection nasale – Mode d'emploi
ES Protección nasal – Instrucciones de uso
IT Protezione nasale – Istruzioni per l'uso
NL Neusbescherming – Gebruiksaanwijzing
GR Προστασία της μύτης – Οδηγίες χρήσης
TR Burun koruması – Kullanım Talimatları
PL Ochrona nosa – Instrukcja obsługi
BG Защита за носа– Инструкции за употреба
CZ Nosní ochrana – Návod k použití
HU Orrvédelem – Használati útmutató
LT Nosies apsauga – Naudojimo instrukcijos
NO Nesebeskyttelse – Bruksanvisning
RO Protecția nasului – Instrucțiuni de utilizare
SK Ochrana nosa – Návod na použitie
SE Nässkydd – Bruksanvisning

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GB CPAPNose Instructions for Use
CPAPNose by KREIENBAUM
Item numbers and available sizes
• KRE.0164S CPAPNose small 04260642890018
• KRE.0165L CPAPNose large 04260642890025
Intended use of the product CPAPNose is a hydrocolloid nasal strip for skin protection against pressure sores and skin lesions during non-invasive CPAP ventilation with a mask.
Special storage and/or handling conditions For optimal storage conditions, the temperature should be in the range from 10 to 25°C. Store in a dry place, protected from direct sunlight. CPAPNose is to be applied to the patient immediately after opening the individual packaging.
Note on sterility and sterilization The product is non-sterile. It is manufactured under clean-room conditions.
Note on single use product CPAPNose is a medical device for single use and can be used on a single patient for up to five days.
Details on indications, contraindications Skin protection of premature babies during CPAP ventilation. Contraindications: Use only permitted on fully developed skin.
Patient target group(s) Premature and newborn babies.
Intended users For use by healthcare professionals or under their supervision (doctor, midwife or nurse) in the operating room, recovery room, anesthesia room, intensive care unit, emergency room and/or in medical, obstetric, maternity and surgical departments.
Expected clinical benefit Prevention of pressure marks in the nasal area during CPAP ventilation of premature babies.
Preparation Clean and degrease the skin. Do not use alcohol. Make sure that the skin is clean and dry. The skin should not have been treated with oils or skin lotions, as these can impair the adhesive effect of CPAPNose.
Application The hydrocolloid must be actively heated before it is removed from the plastic film. Warm the CPAPNose in your hand for approx. 60 seconds, for example, or place it under a radiant warmer for 15 seconds or hold it against a chemical heel warmer for 10-15 seconds. Remove the CPAPNose from the plastic film/base plate. Apply the CPAPNose starting at the base of the nasal bone. Carefully stroke the CPAPNose over the bridge of the nose. Attach the narrow saddle to the nasal septum. Press the wider saddle under the nose (above the lip). Let the palm of your hand rest on the CPAPNose for a few seconds. The heat is required for the product to adhere well to your skin (approx. 60 seconds). Then the CPAP mask / prong may be used as usual. The mask should rest on the CPAPNose.
Using the CPAPNose without a nose saddle
The hydrocolloid must be actively heated before it is removed from the plastic film. Warm the CPAPNose in your hand for approx. 60 seconds, for example, or place it under a radiant warmer for 15 seconds or hold it against a chemical heel warmer for 10-15 seconds. Take a sterile pair of scissors and cut off the nose saddle on the wide side of the CPAPNose. The result is a triangle without a saddle. The removed part may no longer be used and is disposed of in accordance with hospital regulations. Remove the CPAPNose from the plastic film/base plate. Apply the CPAPNose starting at the base of the nasal bone. Carefully stroke the CPAPNose over the bridge of the nose. Let the palm of your hand rest on the CPAPNose for a few seconds. The heat is required for the product to adhere well to the skin (approx. 60 seconds). Then the CPAP mask / prong may be used as usual. The mask should rest on the CPAPNose.
Removal Use lukewarm water or saline solution to remove the CPAPNose. Carefully remove one corner of the CPAPNose with a swab soaked in water and slowly pull the CPAPNose from your nose using the soaked swab. Then dispose of the product according to hospital regulations.
Precautions and warnings After opening the individual packaging, the product must be used immediately to ensure sterility. Remove immediately if skin irritation occurs. Moist skin and skin treated with lotions or oils impair the adhesive effect of the CPAPNose. Do not use alcohol to clean the skin before use.
Serious incidents Any serious incident that has occurred in connection with this medical device should be reported to KREIENBAUM Neoscience GmbH and the competent authority of the Member State in which the user and/or patient is located.

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